



HOW TO COMPLETE AND SUBMIT AED INSPECTION

Go online to <https://secure.heart-safe.ca/default.cfm>


Login to the Member Sign-In with your username and password



Click on the AED Inspection Icon



Click on the AED you want to complete an Inspection for



AED Inspection
home - member area - aed inspection


Please ensure an Inspection is submitted for each AED every 30 days. Inspections may be submitted at any time.
Click on each AED record to complete the Inspection.

AED Inspection Checklist(s)			
SITE	AED NICKNAME/MODEL	SERIAL NUMBER	LAST INSPECTION
testsite1	✓ test1 Zoll AED Plus	GOOG000000	28 days ago
testsite1	🚩 Philips HeartStart OnSite	SSSSSSSSSS	35 days ago

[Click to submit a checklist for this AED](#)

Continued on page 2.....

1. Ensure all boxes are checked off
2. Optional comments for your reference can be entered here, after the box is checked
3. Click submit to log your inspection



AED Inspection
home - member area - aed inspection

AED Nickname	AED Model	Serial Number	Site
...	Philips HeartStart OnSite	SSSSSSSSSS	testsite1

A. AED

- AED Readiness indicator displaying green or "OK". 07-06-2016
- AED pads or cartridge attached in unopened pouch and free of damage or foreign debris. 07-06-2016
- Battery Pack installed in AED. 07-06-2016
- No "Errors" or "Alarms" noted. 07-06-2016
- Spare Pads in AED carrying case pouch or under cover. 07-06-2016
 Click here if not equipped.
- Pediatric Pads/Key/Switch Present. 07-06-2016
 Click here if not equipped.

B. RESPONDER KIT INCLUDES
Kit should include: Pocket mask or shield (1), Disposable gloves (1), Razor (1), Scissors (1), Towel (1)


- Kit Complete and attached to the AED. 07-06-2016

C. ALARM BOX

- Check Cabinet alarm and replace battery every 12 months. 07-06-2016
 Click here if not equipped.

Submit Cancel

You will get confirmation the inspection has been submitted



SUCCESS: YOUR AED INSPECTION HAS BEEN SUBMITTED.

Submit Cancel

Your AED summary shows the inspection submitted today

AED Inspection
home - member area - aed inspection

Please ensure an Inspection is submitted for each AED every 30 days. Inspections may be submitted at any time.
Click on each AED record to complete the Inspection.

AED Inspection Checklist(s)			
SITE	AED NICKNAME/MODEL	SERIAL NUMBER	LAST INSPECTION
testsite1	✓ test1 Zoll AED Plus	GOOG000000	28 days ago
testsite1	✓ Philips HeartStart OnSite	SSSSSSSSSS	today

END